




# TrueBond

# One Centrifuge for all Your Veterinary Needs



Using Drucker Diagnostics' new, custom-designed centrifuge, veterinary laboratories can now spin microhematocrit tubes, fecal, blood, urine and more in one centrifuge.



## TrueBond

Multipurpose Veterinary Centrifuge

- Capacity: 6 x 3-15 mL tubes, 2 x .5-1.5 mL tubes, 8 x microhematocrit tubes
- Speed: 1850 xg, 3900 RPM
- Applications: Microhematocrit, fecal, blood, urine



**Our Promise:** Made in the USA. Fast order fulfillment. Unmatched reliability and performance. Great support.